

FERTILITY



Investigation of  
**FEMALE** REPRODUCTIVE  
HORMONE **DYSFUNCTIONS**

from diagnosis,  
the seeds of better health



# contents

Precocious and delayed puberty

Secondary amenorrhea

Hirsutism

Other pathologies : primary amenorrhea  
and hyperprolactinemia

Menopause

**The approach** used for each of the subjects discussed  
in this booklet is identical :

- brief physiological description
- clinical approach
- basic biological profile
- interpretation of results
- secondary examinations, if required
- treatment


Dynamic tests and a list of the main hormone assays  
are given at the end of the booklet.

# early follicular phase


D-3 ⇨ D5 = selective follicular recruitment

# late follicular phase

# luteal phase





OVULATION



Steroidogenesis : LH-dependent / Aromatization : FSH-dependent


- theca interna
- theca externa
- vascularization
- granulosa cells
- E<sub>2</sub> ..... Estradiol
- Pg ..... Progesterone
- GnRH ..... Gonadotropin Releasing Hormone
- FSH, LH ... Gonadotropins
- IGF ..... Insulin-like Growth Factor
- PGE, PGF ... Prostaglandins



HORMONAL PHYSIOLOGY

# puberty

biochemical mechanisms



Complete pubertal development takes 2 to 3 years. It is preceded by an adrenal maturation phase (at the age of 6 or 7) known as the adrenarche, biochemically characterized by an increase in circulating DHEAS\*.

\*Dehydroepiandrosterone sulfate.

# precocious puberty

Onset of puberty before the age of 8 (European population).



## CLINICAL SIGNS

Breast development and/or growth of pubic and axillary hair.



## INITIAL PROFILE

This profile aims to differentiate between :

- > isolated pubic and axillary hair growth (pubarche)
- > isolated breast development (thelarche)
- > central precocious puberty
- > primary precocious puberty (pseudoprecocious puberty)

It comprises :

- > basic FSH and LH levels + LH-RH test (GnRH)
- > Estradiol
- > DHEAS to evaluate adrenal maturation or adrenarche
- > evaluation of stature and bone age



## INTERPRETATION OF RESULTS

secondary sexual characteristics

basic FSH-LH levels

response to LH-RH test

secondary examinations for confirmation or orientation

dynamic tests

	Isolated or predominant breast development	Predominant or isolated axillary hair growth	More or less balanced development
secondary sexual characteristics	isolated breast development	isolated axillary hair growth	balanced development
basic FSH-LH levels	FSH ↑, LH low	LH ↑, FSH normal or low	normal or increased
response to LH-RH test	FSH ↑, LH prepubertal or low	LH ↑, FSH prepubertal	pubertal
secondary examinations for confirmation or orientation	<p>PRECOCIOUS THELARCHE</p> <p>OVARIAN PRIMARY PRECOCIOUS PUBERTY (OR PSEUDOPRECOCIOUS PUBERTY)</p> <ul style="list-style-type: none"> <li>abdominal-pelvic radio-imaging techniques</li> <li>associated endocrine disorders</li> <li>skeleton radiography</li> </ul>	<p>adrenal disorders</p> <ul style="list-style-type: none"> <li>abdominal-pelvic radio-imaging techniques</li> <li>4-Androstenedione, testosterone, DHEAS, 17-OH Progesterone</li> </ul>	<p>TRUE PRECOCIOUS PUBERTY (OR CENTRAL)</p> <ul style="list-style-type: none"> <li>cerebral radio-imaging techniques</li> </ul>
dynamic tests	<p>ovarian tumor</p> <p>functional cysts</p> <p>McCune-Albright syndrome</p> <p>PCO</p>	<p>PRECOCIOUS PUBARCHE (ADRENARCHE)</p> <p>ADRENAL GLAND PRIMARY PRECOCIOUS PUBERTY (OR PSEUDOPRECOCIOUS PUBERTY)</p> <p>synacthene, dexamethasone</p> <p>tumor</p> <p>Cushing's syndrome</p> <p>hyperplasia (21-hydroxylase deficiency)</p>	<p>CENTRAL NEUROGENIC OR IDIOPATHIC PRECOCIOUS PUBERTY</p>



## TREATMENT

- In cases of true central precocious puberty, pubertal development is halted using an LH-RH agonist (an annual LH-RH test controls the degree of pituitary blockage).

- Treatment of congenital adrenal hyperplasia (CAH).
- Treatment of the tumor, if required.

# delayed puberty

No signs of puberty after the age of 13 - 14 (European population).



## CLINICAL SIGNS

No breast development, nor pubic and axillary hair growth.



## INITIAL PROFILE

This profile aims to differentiate between :

- > delayed puberty
- > hypogonadotropic hypogonadism
- > hypergonadotropic hypogonadism

It comprizes :

- > basic FSH and LH levels + LH-RH test (GnRH)
- > Estradiol
- > DHEAS to evaluate adrenal maturation or adrenarche
- > evaluation of stature and bone age for orientation of diagnosis to delayed puberty
- > assay of prolactin to eliminate hyperprolactinemia




## INTERPRETATION OF RESULTS

secondary sexual characteristics

basic FSH-LH levels

response to LH-RH test

secondary examinations for confirmation or orientation



## TREATMENT

substitutive (estrogen then estrogen-progestrone), except in cases of delayed puberty.

# secondary amenorrhea

(without hirsutism) Women < 46 yrs old



## CLINICAL SIGNS

- no specific clinical signs
- no menses for over 3 months

### Anamnesis :

Date of last childbirth, variation in weight, drugs, genital and breast examination, stop estrogen-progesterone treatment, affective shock...

Background history of :

- chemotherapy and radiotherapy.
- surgery (ectopic pregnancy, ovariectomy, appendicectomy...).
- infection (salpingitis, STD, tuberculosis)



## INITIAL PROFILE

### Firstly :

- > assay hCG to exclude pregnancy
- if hCG negative, make an appointment in 1 to 2 weeks time using a menothermal curve (to exclude pregnancy or trophoblastic tumor)

### Then :


- > FSH, Estradiol, (LH)
  - > Prolactin (PRL)
  - > TSH if apathy and/or weight gain.
- Hypothyroidism (increased TSH) leads to an increase in TRH which stimulates PRL secretion.



## INTERPRETATION OF RESULTS

hCG : negative  
TSH : normal or regulated

### Basic tests



### secondary examinations for confirmation and orientation

- iatrogenous origin
- adenoma (radio-imaging techniques of the sella turcica and visual field)

Additional investigations and/or action to be taken :  
**Either desire for pregnancy :**  
- FSH, LH, E<sub>2</sub>  
- Ovulation induction  
**Or wait :**  
- For menses to return  
- FSH and LH control  
- Cyclic progesterone treatment, if required

**Etiologies**  
- Infection  
- Endometrium lesion (following curettage)  
- Physical / Psychological trauma (central amenorrhea)  
- Sport at competitive level  
- Adenoma / adrenal tumor or secreting ovary ? E<sub>2</sub>, Cortisol, Testo  
- Recent non-secreting pituitary adenoma  
- Hemorrhagic childbirth (Sheehan's syndrome)

**Depending on etiology**  
- Hysteroscopy  
- Laparoscopy  
- Endometrium, cervix biopsy  
- Pelvic radio-imaging techniques  
- Scan, NMR

# Secondary amenorrhea / spaniomenorrhea with hirsutism



## CLINICAL SIGNS

No menses for over 3 months or 2 to 4 cycles per year.

- hirsutism : excess hair growth in regions stimulated by sexual hormones.  
Possible to grade (0 to 4) the level of excess hair.  
Hirsutism is pathological while hypertrichosis is ethnic and family-related
- acne, seborrhea
- possible obesity (android fat distribution ?)
- recent signs of virilization (voice deepening, clitoromegaly...)

### Anamnesis :

- date of puberty - regularity of menses
- weight gain
- evolution of hirsutism - treatment in progress
- menothermic curves



## INITIAL PROFILE (before D5 if spaniomenorrhea)

This profile aims to distinguish :

- > the origin of hyperandrogenia (ovarian, adrenal or idiopathic).

Cases of recent virilization may also be of tumor origin, for which biological diagnosis must be rapid.

It comprizes :

- > FSH, LH, Prolactin (PRL),
- > Testosterone (T), E<sub>2</sub>, 4-Androstenedione (Δ4), DHEAS, 17-OH-Progesterone (17-OH-P)
- > Free cortisol in urine



## INTERPRETATION OF RESULTS

FSH	N	N	N	N	N
LH	↗↗	N or ↗	N or ↗	N or ↗	N
LH/FSH ratio (before D5) N<1	>2	N	N or >1	non informative	non informative
4-Androstenedione	↗↗	N or ↗	N or ↗	↗	↗↗↗
Testosterone	N or ↗	N or ↗	N greater or ↗	N	↗↗
		bio-available T ↗			
DHEAS	N or ↗	N	↗	↗	N
free cortisol in urine/creatinine ratio	N	N	N	↗↗	N
PRL	N or ↗	N	N	N	N
E <sub>2</sub>	↗	N	non informative	non informative	↗
17-OH-P	N	N	N greater or ↗	↗	N or ↗
	<b>PCO</b> Polycystic ovarian syndrome	Idiopathic hirsutism, hyperproduction and hyperconsumption of androgens, 5-α-reductase-receptor disease	Late appearance of adrenal hyperplasia (21-hydroxylase deficiency)	Specific investigation of adrenocortical hyperfunction (Cushing's syndrome)	Virilizing ovarian tumor
<b>Dynamic tests</b>					
LH-RH test	LH ↗↗ FSH normal	N	N		
ACTH test	-	N	17-OH-P > 5 ng/ml at T <sub>0</sub> + 60 mins		-
Radio-imaging techniques	Ovarian cysts	Normal ovaries	Adrenal volume ↗		- scan, NMR - ovarian vein sampling : assay of T and Δ4 for each ovary

N = reference value

# primary amenorrhea

adolescents, young women



## CLINICAL SIGNS

Secondary sexual characteristics are most often absent (impuberty)

- external genital organs (pubic hair, vulva, clitoris) not or only slightly developed

- internal genital organs (vagina, uterus, hymen, ovaries) normal or more or less absent



## INITIAL PROFILE

FSH, LH, E<sub>2</sub>



## INTERPRETATION OF RESULTS

### Basic tests

FSH ↑↑, LH ↑↑, E<sub>2</sub> ↓

**HYPERGONADOTROPIC  
HYPOGONADISM**

FSH, LH = N or ↓  
E<sub>2</sub> low

**HYPOGONADOTROPIC  
HYPOGONADISM**

**DELAYED  
PUBERTY**

### secondary examinations for confirmation or orientation

- sex chromatin (jugal mucous)
- karyotype

- LH-RH test

**46XY**  
<10% of cases

**Swyer-James  
syndrome**

- no breasts
- no testicles
- streak ovaries

**Testicular Feminizing  
syndrome**

- FSH ↑ (or N)  
LH ↑↑  
Testo type ♂
- Minimum breast development
  - Abdominal or inguinal testes = androgen insensitivity (no DHT-receptors)

**46X0**  
90% of cases

**Turner's  
syndrome**

- If LH-RH test ovaries in strip + nanism more or less variable malformations

**46XX**  
<10% of cases

- Gonadal dysgenesis
- Gynatresia, breast and hair growth ± N

**17-α-hydroxylase  
deficiency**

- no cortico-adrenal and sexual steroidogenesis
- + High blood pressure
- + alkalosis

+

**hypothalamic  
origin**

- Kallman's syndrome
- Prader-Willi's syndrome
- Functional : coeliac disease, chronic disease (IRC...) anorexia

-

**pituitary origin**

- infiltration disease, sarcoidosis, histiocytosis
- Tumors : cranial-pharyngioma glioma adenoma



# hyperprolactinemia



## PHYSIOLOGICAL ROLE - REGULATION

Prolactin (= lactotropic hormone) is a peptide hormone

**Role** : triggering and maintaining lactation

**Secretion** by the pituitary according to a circadian rhythm

(max. : between 2 am and 6 am ; min. : 10 am and 12 pm)

### Regulation :

*inhibition* of secretion by hypothalamic Dopamin

*stimulation* of secretion by TRH (see dynamic tests page 20)

other factors increase secretion : stress, thyroid hormones, corticoids, estrogens, physical exercise, meals and some therapeutic drugs (see below).



## INITIAL PROFILE

### Precautions when assaying Prolactin :

*rest* (no stress)

*detailed questioning* : dates of last menses, intensive sport, treatment in progress...

### Indicative reference values :

men and children : < 15 ng/ml

menstruating women : < 25 ng/ml

menopausal women : < 20 ng/ml

pregnant women : progressive increase until childbirth



## INTERPRETATION OF RESULTS

### Hypoprolactinemia

There is no hyosecretion threshold ; low levels have no clinical significance.

### Hyperprolactinemia

#### clinical signs

#### basic prolactin

response to TRH test

#### secondary examinations

#### diagnosis

	primary	secondary	functional	biological
	cephalea, amenorrhea, galactorrhea		infertility, menstruating disorders	no clinical or radiological signs
	> 120 ng/ml	50 - 100 ng/ml	approx. 50 ng/ml	> 30 ng/ml
	no response			
	CT Scan, NMR			Chromatography of different circulating PRL forms
	pituitary micro or macroadenoma	hypothalamus control suppressed due to tumor lesion, or post radiotherapy or sarcoidosis	following endocrine pathologies : hypothyroidism, PCO....	secretion of a minor form of PRL : <b>BIG-BIG Prolactin</b> ; no pathology



## TREATMENT

In cases of primary hyperprolactinemia, treatment using Dopamin agonist and / or surgery.

### Therapeutic drugs causing hyperprolactinemia (examples)

**Psychotropic drugs**

**Antidepressant drugs**

**Estrogens**

**Metoclopramide**

**Opium-based drugs**

**Cimetidine**

**Romitidine**



### CLINICAL SIGNS

**Perimenopause** : physiological situation as of 45 - 50 yrs of age, with irregular or shorter cycles: follicular phase increasingly shorter (FSH  $\uparrow$ ,  $E_2$   $\downarrow$ ), then luteal insufficiency with Pg  $\downarrow$  and LH  $\uparrow$  which stimulates  $E_2$ .  
Accompanying clinical signs : mastodynia, abdominal and pelvic distension, weight gain, irritability, hot flushes.

**Confirmed menopause** : estrogen levels markedly reduced (no menses for over one year).

**Early menopause** : after surgical or chemical castration, irradiation, chemotherapy, intense stress, affective shock, pituitary adenoma, pituitary neurosurgery...



### INITIAL PROFILE

**Perimenopause** : FSH before D5  $\uparrow$   
(decrease of follicular stock).

**Confirmed menopause** :  $E_2$   $\downarrow\downarrow$ , FSH  $\uparrow\uparrow$ , LH  $\uparrow$



### TREATMENT AND/OR FOLLOW-UP

The strong hormonal instability (both intra- and inter-individual) most often requires progesterone treatment until menstruation stops.

#### Substitutive hormone treatment

Short- and long-term beneficial effects :

- Quality of life
- Prevention of heart and cardio-vascular disease
- Prevention from the risk of osteoporosis
- Evaluation of the benefits to the heart and bones *versus* the risk of breast cancer

#### Biological monitoring of substitutive hormone treatment

- If required,  $E_2$  to adapt the posology if percutaneous substitutive hormone treatment
  - FSH (<30 IU/l)
  - Carbohydrate-lipid profile annually or every 2 yrs.
- Essentially clinical and radiological monitoring (mammography, ultra-sound scan).


# dynamic investigation tests




## LH-RH (or GnRH) TEST

Intravenous **injection** at  $T_0$ , of 100  $\mu\text{g}/\text{m}^2$  (child) or 100  $\mu\text{g}$  (adult) of LH-RH

**Assay** of FSH or LH at  $T_{-15}$ ,  $T_0$ ,  $T_{20}$ ,  $T_{40}$ ,  $T_{60}$ ,  $T_{90}$  minutes



prepubertal type response



pubertal type response



## L-DOPA OR TRH TEST

Production of Prolactin is

- reduced by L-dopa

Assay of PRL at  $T_0$ ,  $T_{15}$ ,  $T_{30}$ ,  $T_{60}$ ,  $T_{90}$  and  $T_{120}$  minutes  
maximum slowing down at  $T_{60}$  or  $T_{90}$

- stimulated by TRH

Assay of PRL at  $T_0$ ,  $T_{15}$ ,  $T_{30}$ ,  $T_{60}$ ,  $T_{90}$  and/or  $T_{120}$  minutes  
increase of 200 to 300% between  $T_{15}$  and  $T_{60}$



## SYNACTHEN TEST (SYNTHETIC ACTH)

**Intramuscular injection** of 0.25 mg of Synacthen at  $T_0$  (0.125 mg if  $\leq 2$  yrs old).

**Assay** of cortisol, 17-OH-progesterone, aldosterone, DHEAS, 4-Androstenedione at  $T_0$ ,  $T_{30}$  and/or  $T_{60}$  minutes.

An objective normal response is obtained if there is an increase in cortisol and aldosterone (minimum factor 2), without any significant modification of the other parameters.

## **BLOOD HORMONE ASSAYS :**

---

<b>VIDAS hCG</b> .....	<b>ref. 30 405</b>
<b>VIDAS LH</b> .....	<b>ref. 30 406</b>
<b>VIDAS FSH</b> .....	<b>ref. 30 407</b>
<b>VIDAS Prolactin</b> .....	<b>ref. 30 410</b>
<b>VIDAS Progesterone</b> .....	<b>ref. 30 409</b>
<b>VIDAS Estradiol II</b> .....	<b>ref. 30 431</b>
<b>VIDAS Testosterone</b> .....	<b>ref. 30 418</b>
<b>VIDAS Cortisol</b> .....	<b>ref. 30 417</b>
<b>VIDIA hCG*</b> .....	<b>ref. 38 300</b>
<b>VIDIA LH*</b> .....	<b>ref. 38 310</b>
<b>VIDIA FSH*</b> .....	<b>ref. 38 320</b>
<b>VIDIA Prolactin*</b> .....	<b>ref. 38 330</b>
<b>VIDIA Progesterone*</b> .....	<b>ref. 38 340</b>
<b>VIDIA Estradiol*</b> .....	<b>ref. 38 350</b>

Availability of some VIDAS tests may be restricted in certain countries due to registration requirements. Consult our local representatives for further information.

\*In development.

This booklet was created in collaboration with :

**Dr. I. COLLIGNON**  
(Versailles Hospital Center, FRANCE)

**Prof. D. PORQUET**  
(R. Debré University Hospital Center, Paris, FRANCE)

**bioMérieux sa**  
69280 Marcy l'Etoile  
France  
Tel. : 33 (0)4 78 87 20 00  
Fax : 33 (0)4 78 87 20 90  
[www.biomerieux.com](http://www.biomerieux.com)

