

FERTILITY



Investigation of
MALE REPRODUCTIVE
HORMONE **DYSFUNCTIONS**

from diagnosis,
the seeds of better health



contents

Precocious and delayed puberty

Hypogonadism

Gynecomastia

Azoospermia



The approach used for each of the subjects discussed in this booklet is identical :

- brief physiological description
- clinical approach
- basic biological profile
- interpretation of results
- secondary tests, if required
- treatment

Dynamic tests and a list of the main hormone assays are given at the end of the booklet.

gonadotropic axis (adult)

gonadotropic axis (puberty)



Complete pubertal development takes 2 to 3 years. It is preceded by an adrenal maturation phase (at the age of 7 or 8), known as the adrenarche, biochemically characterized by an increase in circulating DHEAS*.

* Dehydroepiandrosterone sulfate.

precocious puberty

Onset of puberty before the age of 9
(European population)



CLINICAL SIGNS

Growth of pubic hair, and possible signs of gynecomastia.

Growth of testes and/or development of the penis and scrotum.



INITIAL PROFILE

This profile aims to differentiate between :

- > central precocious puberty.
- > primary precocious puberty or pseudoprecocious puberty.

It comprizes :

- > basic FSH and LH levels + LH-RH test (GnRH).
- > Testosterone.
- > DHEAS to evaluate adrenal maturation or adrenarche.
- > evaluation of stature and bone age.



INTERPRETATION OF RESULTS

| | 2 small testes | Only 1 enlarged testis | 2 enlarged testes |
|--|---|---|---|
| secondary sexual characteristics | | | |
| basic FSH-LH levels | <i>low</i> | <i>low</i> | <i>normal or increased</i> |
| response to LH-RH test | <i>prepubertal</i> | <i>prepubertal or low</i> | <i>pubertal</i> |
| | ADRENAL PRIMARY PRECOCIOUS PUBERTY (or pseudoprecocious puberty) | TESTICULAR PRIMARY PRECOCIOUS PUBERTY (or pseudoprecocious puberty) | TRUE PRECOCIOUS PUBERTY (or central precocious puberty) |
| secondary tests for confirmation or orientation | <ul style="list-style-type: none"> • abdominal / pelvic radio-imaging techniques • normal testosterone or \uparrow, DHEAS and $\Delta 4$, 17-OH-P \uparrow | <ul style="list-style-type: none"> • high testosterone level • radio-imaging techniques of the testes | <ul style="list-style-type: none"> • brain radio-imaging techniques • "pubertal" testosterone |
| dynamic tests | <p>ACTH test</p> <ul style="list-style-type: none"> Isolated pubic hair or precocious adrenarche Congenital adrenal hyperplasia Adrenal tumor | Tumor of the testes | Central, idiopathic or neurogenic precocious puberty |



TREATMENT

Treatment of the tumor, if required.
In cases of true central precocious puberty, pubertal development is halted by using an LH-RH agonist (an annual LH-RH test controls the degree of pituitary blockage).

Treatment of congenital adrenal hyperplasia.

delayed puberty

No signs of puberty after the age of 15
(European population)



CLINICAL SIGNS

No growth of pubic hair, nor development of other secondary sexual characteristics.



INITIAL PROFILE

This profile aims to differentiate between :


- > delayed puberty: hCG test.
- > hypogonadotropic hypogonadism.
- > hypergonadotropic hypogonadism.

It comprizes :

- > an hCG test.
- > basic FSH and LH levels + LH-RH test.
- > Testosterone.
- > DHEAS to evaluate adrenal maturation or adrenarche.
- > evaluation of stature and bone age to orient a diagnosis of delayed puberty.



INTERPRETATION OF RESULTS



TREATMENT

Substitutive (testosterone-based) for functional hypogonadotropic hypogonadism and some types of hypergonadotropic hypogonadism.

acquired hypogonadotropic hypogonadism



CLINICAL SIGNS

- decreased libido
- sexual indifference
- impotence
- physical fatigue
- reduced muscle mass
- reduced testicular volume



INITIAL PROFILE

Testosterone, FSH, LH, Prolactin, Inhibin B



INTERPRETATION OF RESULTS

testosterone $\searrow\searrow$ to $\searrow\searrow\searrow$
(for info : Reference values = 3 to 12 ng/ml)

inhibin B \searrow to $\searrow\searrow$
FSH and LH are normal to $\searrow\searrow$

ETIOLOGIES

SUPRASELLAR LESIONS

- Tumoral (CT scan + NMR)
 - craniopharyngioma, visual disorders and/or panhypopituitarism
 - germinoma, (hCG, AFP)
 - glioma of optic chiasm or hypothalamus
- Post-infectious
 - tuberculous meningitis or other microorganism
- Infiltration processes (NMR necessary)
 - histiocytosis
 - sarcoidosis
- Post radiotherapy
- Section of the pituitary stalk

PITUITARY ORIGIN

- Hemochromatosis
 - ferritin > 1000 $\mu\text{g/l}$
 - investigation of C282Y mutation
 - ▶ gonadotropic cells affected by Fe deposition
 - ▶ no GnRH stimulation
- Hyperprolactinemia
 - Prolactin = strong antagonistic effect on the hypothalamus
 - ▶ Prolactin adenoma (80% of pituitary adenomas)
- Secreting or non-secreting pituitary adenomas

FUNCTIONAL DEFICIENCIES

- Anorexia (rare in males)
- Nutritional disorders : celiac disease due to gluten intolerance
- Hypercorticism, long-term corticosteroid therapy
- Hyperestrogenemia
 - estrogen treatment (transvestites)
 - estradiol-secreting tumor

- radio-imaging techniques, GnRH test

secondary tests for confirmation or orientation



TREATMENT

- Pulsatile GnRH pump (1 to 2 years), if positive test response
- Long-term IM administration of testosterone
- Specific according to etiology

gynecomastia



CLINICAL SIGNS

Unilateral or bilateral enlargement of the male breast (normally undeveloped)

- Newborn (due to placental estrogen)
- Puberty (60 % of adolescents aged between 12 and 15)
- Adult : frequent, often asymptomatic, increasing incidence with age, occasionally with breast deformation and tenderness.

Anamnesis :

- hyperthyroidism
- renal insufficiency (dialysis)
- hepatic insufficiency
- therapeutic drugs or non-therapeutic substances *
- libido
- galactorrhea



INITIAL PROFILE

Biological :

Testosterone, Estradiol, hCG

Radiological :

mammography

hCG secretion may orient diagnosis towards a secreting chorionic carcinoma or neoplasia (e.g. lung). Mammography should identify an increase in adipose tissue linked to obesity or breast cancer (rare).



INTERPRETATION OF RESULTS

| | | | | | |
|--------------|---|---|---|---|---|
| Testosterone | ↗ | ↘ | ↘ | N | N |
| Estradiol | N | N | ↘ | N | ↗ |
| FSH-LH | N | ↗ | ↘ | N | - |
| SHBG | - | - | - | ↗ | - |
| Prolactin | N | - | ↗ | - | - |

| | | | |
|--|---|--|---|
| <p>ANDROGEN RESISTANCE</p> <p>▼</p> <p>Testosterone-receptor deficiency</p> <p>▼</p> <p>Reifenstein's syndrome</p> | <p>PERIPHERAL HYPOGONADISM</p> <p>Primitive karyotype (Klinefelter's syndrome)</p> <p>Secondary</p> <ul style="list-style-type: none"> - viral orchitis - castration - sarcoidosis | <p>HYPOGONADOTROPIC HYPOGONADISM</p> <p>Pituitary adenoma</p> <ul style="list-style-type: none"> - with PRL secretion (galactorrhea) - other | <p>SECRETING ESTRADIOL TUMOR</p> <p>• Ultrasonography of the testes</p> <p>Leydig cell tumor</p> <p>• Adrenal CT scan</p> |
|--|---|--|---|

N = reference value

* Therapeutic drugs responsible for gynecomastia

- Estrogens
- Androgens, anabolic steroids (peripheral aromatizations)
- Spironolactone
- Cimetidine
- hCG therapy
- Antiandrogenic drugs
- Digitalis
- Isoniazid

- Cytostatics and irradiation (testicular lesion)
- Hyperprolactinemia drugs
- Neuroleptics, tricyclic antidepressant drugs

* Non-therapeutic substances

- Cannabis (marijuana)
- Heroin
- Estrogen-rich hair lotions
- Partner's local estrogen therapy (« vacation gynecomastia »)

azoospermia

The diagnosis of azoospermia is based on the absence of spermatozoa.



CLINICAL SIGNS

The clinical approach takes into consideration :

- the size and consistency of the testes
- the possible existence of gynecomastia
- hypoandrimism
- mental retardation
- pubic and axillary hair and muscle mass
- size of penis (micropenis)




INITIAL PROFILE

Involves FSH, LH and Testosterone



INTERPRETATION OF RESULTS



secondary tests for confirmation or orientation

Reminder : characteristics of a normal spermogram

(histological viewpoint)


- Ejaculate volume > 2.5 ml
- Number of spermatozoa : > 20 million/ml
- Percentage of motile spermatozoa after 4 hrs > 60 %
- Percentage of morphologically normal forms 40 - 70 %
- Percentage of living forms (vitality) 60 - 90 %




LH-RH (or GnRH) TEST

Intravenous injection at T_0 of 100 $\mu\text{g}/\text{m}^2$ (child) or 100 μg (adult) of LH-RH

Assay of FSH or LH at T_{-15} , T_0 , T_{20} , T_{40} , T_{60} , T_{90} minutes



prepubertal type response



pubertal type response



hCG TEST

Intramuscular injection on D_1 , D_3 , and D_5 of 1500 IU of hCG,

Testosterone assay on D_1 and D_6 .

The testosterone level must be at least 3 ng/ml in prepubertal phase and may reach similar values to those of the adult (7 to 12 ng/ml) during puberty.



SYNACTHEN TEST (synthetic ACTH)

Intramuscular injection at T_0 of 0.25 mg of Synacthen (0.125 mg if ≤ 2 years old),

Assay of cortisol, 17-OH-progesterone, aldosterone, DHEAS, 4-Androstenedione at T_0 , T_{30} and/or T_{60} minutes.

An objective normal response is obtained if there is an increase in cortisol and aldosterone (minimum factor 2) without any significant modification of the other parameters.

BLOOD HORMONE ASSAYS :

| | |
|----------------------------------|--------------------|
| VIDAS hCG | ref. 30 405 |
| VIDAS LH | ref. 30 406 |
| VIDAS FSH | ref. 30 407 |
| VIDAS Prolactin | ref. 30 410 |
| VIDAS Progesterone | ref. 30 409 |
| VIDAS Estradiol II | ref. 30 431 |
| VIDAS Testosterone | ref. 30 418 |
| VIDAS Cortisol | ref. 30 417 |
| | |
| VIDIA hCG* | ref. 38 300 |
| VIDIA LH* | ref. 38 310 |
| VIDIA FSH* | ref. 38 320 |
| VIDIA Prolactin* | ref. 38 330 |
| VIDIA Progesterone* | ref. 38 340 |
| VIDIA Estradiol* | ref. 38 350 |

Availability of some VIDAS tests may be restricted in certain countries due to registration requirements. Consult our local representatives for further information.

*In development.

This booklet was created in collaboration with :

Dr. I. COLLIGNON
(Versailles Hospital Center, FRANCE)

Prof. D. PORQUET
(R. Debré University Hospital Center, Paris, FRANCE)

bioMérieux sa
69280 Marcy l'Etoile
France
Tel. : 33 (0)4 78 87 20 00
Fax : 33 (0)4 78 87 20 90
www.biomerieux.com

